

闪 胶

See Flash Run, More
5 minute DNA separation.

Separate DNA, and now RNA, faster than ever with **FlashGel System**

- Get results in less than 5 minutes
- Watch DNA migrate in real-time - without UV light
- Separate DNA – 10 bp to 4 kb
- Separate RNA – 0.5 kb to 9 kb
- Run up to 34 samples on a single cassette

See Demo Video >
QuickTime and Windows Media formats

闪胶系统(The FlashGel System)

美国 Lonza 公司的 FlashGel 产品线，包括 5 分钟出结果的 DNA 电泳预制胶（有多种型号）和 配套的迷你电泳仪，FlashGel 上样染液和 DNA/RNA 标准等。这是迄今为止世界上最快速的电泳装备。痕量快速 DNA 检测电泳胶，是许多实验室必备的工具。

- 1) DNA 电泳 5 分钟出结果，快速没有第二家可比，可以不时用于急需
- 2) 痕量，有时因为酶切量少且片段短，普通琼脂糖胶难以检测到，就可以用这种特殊预制胶
- 3) 方便省事，不需自己花时间配置凝胶，不需要自己配置电泳缓冲液，不需要用自己的电泳槽【FlashGel:可在五分钟内测出 DNA 及片断】

有关产品信息如下: **electrophoresis**

The FlashGel System is the fastest way to separate DNA and the only way to watch DNA migration as it happens. This revolutionary new tool separates DNA in 2-7 minutes. Monitor gel runs right at the bench, without UV light. The highly sensitive system allows a 5X reduction in DNA levels? saving both time and money.

- * Get results in 5 minutes or less.
- * Watch DNA migrate at your bench, in real-time without UV light.
- * Enjoy outstanding resolution and exquisite sensitivity.
- * See NEW formats for small fragment DNA and high sample

throughput. The FlashGEL System consists of enclosed, disposable, precast agarose gel cassettes and a combination electrophoresis and transilluminator unit.

* FlashGel Cassettes contain precast, prestained agarose gels and buffer?no need for gel preparation, buffer addition or gel staining.

* The FlashGel Dock is an electrophoresis apparatus with a built-in transilluminator that provides both separation and detection.

产品订购信息如下:

57026	FlashGel DNA Starter Kit	Kit
包含	1)57025 FlashGel Dock,EA 2)57023/57031 FlashGel DNA Cassette, 9/pk 3)50472 FlashGel DNA Marker(100-4000bp),150ul 4)50463 FlashGel Loading Dye (sample size),1ml	
57025	FlashGel® Dock	Kit
50473	FlashGel DNA Marker,100-4kb,	500ul
57033	FlashGel DNA Marker 50bp-1.5kb,	500ul
57034	FlashGel DNA Marker 100bp-3kb,	500ul
50462	FlashGel Loading Dye 5xconcentrate,	5x1mL
57023	FlashGel DNA Cassette 1.2%,12+1well,single tier	Kit
57031	FlashGel DNA Cassette 2.2%,12+1well,single tier	Kit

北京毕特博生物技术有限责任公司

网 址 : <http://www.bitebo.com>

电话: 010-82015225

服务热线: 400-833-9299

邮箱: info@bitebo.com

